

# SR Scales®

*“The Measure of Excellence”*

Avian Scale  
Model **SR**®V315

Whenever it is time to look to the sky, let the SRV315 not pass you by. This scale permits veterinarian and zoological specialists to weigh birds up to 3,000 grams or 6 pounds.

The large base, along with its 5 inch perch, permits the bird to be comfortably placed for proper weighing. The bottom of the base also has a bracket that attaches the scale to a table or countertop.

The perch attaches to the bottom of the base so the bird will not be able to move the perch around.

The SRV315 can also function as a gram scale. Simply remove the perch for weighing smaller animals. This scale may also be used as a daily nutritional scale.

Whatever your weighing needs may be, **SR Scales**® is *“The measure of Excellence.”*

**Display Resolution-** 3,000g X 1g

**Display Units-** Grams or Ounces

**5 Inch Perch -** Allows tails to be elevated off the ground

**Battery Operated -** Two (2) CR2032 Lithium Batteries

**Hold-** Freezes the weight on the display

**Auto-Power Shutdown-** Approximately Two (2) Minutes




## SYSTEM SPECIFICATIONS

Max Weight	3000g / 105oz
Display Resolution	1g
Accuracy	3,000g X 1g
Display Units	Grams/Oz
Removable Perch	Perch screws in the bottom of the plate to allow use as a small animal scale or dietary portion scale
Easy to Clean	The materials of the perch will not cause fluid back up
Auto Power Down	About Two (2) Minutes
Bracket	Scale can be attached to table or hung on a wall
Power	Two (2) CR2032 Lithium Batteries
Hold	Will Recall Last Stored Weight
Warranty	Two (2) Years

# Avian Scale

## Model **SR**®V315

### DIMENSIONS in (cm)



### PACKAGING

